

MANUFACTURED BY

**KAY ESS TRADERS**

Email: [ravin.kayestraders@gmail.com](mailto:ravin.kayestraders@gmail.com)


382-A, Gali No.4, Bagh Kare Khan, Kishan Ganj, Delhi-7

Mob.: 9811256040

**INDOOR AED CABINET, MODEL: FL**

Kay Ess AED cabinets feature combination audible and Strobe light alarms. They are made of sturdy heavy-gauge steel. AED Defibrillator Cabinet With Strobe Light Alarm can proudly display your AED for all to see and have the peace of mind knowing that it is protected from dirt, dust and possibly theft. In addition, the AED wall cabinets help the rescuers locate the AED device and provide an easy way to access it in case of an emergency.

### *Dimension*



MANUFACTURED BY


**KAY ESS TRADERS**

Email: ravin.kayestraders@gmail.com

382-A, Gali No.4, Bagh Kare Khan, Kishan Ganj, Delhi-7

Mob.: 9811256040

**INDOOR AED CABINET, MODEL: FL**



**Key features:**

- AED Defibrillator Cabinet With Strobe light Alarm.
- Powder coated steel body for strength.
- The loud alarm emits a sound when the door is opened.
- Clear PC glass window and red & black lettering to facilitate visual identification of the AED.
- Economical Design with Smooth Finish.
- Easy to Install – Four holes are punched in the back of this cabinet for a wall mount.
- Three Side AED Graphics.
- Accommodates most AED's on the market.
- Longest Battery backup; Batteries Included.
- Accessories: Two locking keys along with screws & plastic sleeves for fitting.
- Safety Belt.
- One Key covered with Glass for Emergency
- Hammer for breaking the glass in emergency
- Product Dimension in Inches L x W x H: (6.5x15.5x18)
- Packing Dimension in Inches L x W x H: (9.5x18.5x20.5)
- Colour: White

