

MANUFACTURED BY

**KAY ESS TRADERS**

Email: [ravin.kayestraders@gmail.com](mailto:ravin.kayestraders@gmail.com)

382-A, Gali No.4, Bagh Kare Khan, Kishan Ganj, Delhi-7


Mob.: 9811256040

**INDOOR AED CABINET, MODEL: STAR 5**

**Kay Ess premium AED cabinets feature combination audible and flashing light alarm.**

**They are made of sturdy heavy-gauge steel and are large enough to accommodate additional medical supplies. AED Defibrillator Cabinet with Flashing Light Alarm can proudly display your AED for all to see and have the peace of mind knowing that it is protected from dirt, dust and possibly theft. In addition, the AED wall cabinets helps the rescuers locating the AED device and provide an easy way to access it in case of an emergency.**

### ***Dimension***



MANUFACTURED BY


**KAY ESS TRADERS**

Email: [ravin.kayestraders@gmail.com](mailto:ravin.kayestraders@gmail.com)

382-A, Gali No.4, Bagh Kare Khan, Kishan Ganj, Delhi-7

Mob.: 9811256040

**INDOOR AED CABINET, MODEL: STAR 5**



MANUFACTURED BY

**KAY ESS TRADERS**

Email: ravin.kayestraders@gmail.com

382-A, Gali No.4, Bagh Kare Khan, Kishan Ganj, Delhi-7

Mob.: 9811256040

**INDOOR AED CABINET, MODEL: STAR 5**



**Key features:**

- AED Defibrillator Cabinet with Flashing light Alarm.
- Powder coated steel body for strength.
- The loud alarm emits a sound when the door is opened.
- Clear PC glass window and red & black lettering to facilitate visual identification of the AED.
- Economical Design with Smooth Finish.
- Easy to Install – Four holes are punched in the back of this cabinet for a wall mount.
- Three Side AED Graphics.
- Accommodates most AED's on the market.
- Longest Battery backup; Batteries Included.
- Accessories: Two locking keys along with screws & plastic sleeves for fitting.
- Safety Belt.
- First Aid and compartment in inches L x W x H: (6.5x16.5x8).
- AED Pad compartment.
- Medicine compartment.
- Product Dimension in Inches L x W x H: (6.5x16.5x24.5)
- Packing Dimension in Inches L x W x H: (9.5x19x27.5)
- Colour: White