

MANUFACTURED BY

**KAY ESS TRADERS**

Email: [ravin.kayestraders@gmail.com](mailto:ravin.kayestraders@gmail.com)


382-A, Gali No.4, Bagh Kare Khan, Kishan Ganj, Delhi-7

Mob.: 9811256040

## INDOOR AED CABINET, MODEL: With Lock Without Alarm

Kay Ess AED cabinets without Alarm are made of sturdy heavy-gauge steel. AED Defibrillator Cabinet Without Alarm can proudly display your AED for all to see and have the peace of mind knowing that it is protected from dirt, dust and possibly theft. In addition, the AED wall cabinets helps the rescuers locating the AED device and provide an easy way to access it in case of an emergency.

### *Dimension*



MANUFACTURED BY

**KAY ESS TRADERS**

Email: [ravin.kayess traders@gmail.com](mailto:ravin.kayess traders@gmail.com)

382-A, Gali No.4, Bagh Kare Khan, Kishan Ganj, Delhi-7

Mob.: 9811256040

## INDOOR AED CABINET, MODEL: With Lock Without Alarm

Clear PC Glass

Door  
Handle

Lock

Safety Belt

Key & Hanging  
Accessories

### Key features:

- AED Defibrillator Cabinet Without Alarm.
- Powder coated steel body for strength.
- ClearPC glass window and red & black lettering to facilitate visual identification of the AED.
- Economical Design with Smooth Finish.
- Easy to Install – Four holes are punched in the back of this cabinet for a wall mount.
- Three Side AED Graphics.
- Accommodates most AED's on the market.
- Accessories: Two locking keys along with screws & plastic Sleeves for fitting.
- Safety Belt.
- Lock with Two Keys.
- Product Dimension in Inches L x W x H: (6.5x15.5x16.5)
- Packing Dimension in Inches L x W x H: (9.5x18.5x18.75)
- Colour: White

