

MANUFACTURED BY

KAY ESS TRADERS

Email: ravin.kayestraders@gmail.com


382-A, Gali No.4, Bagh Kare Khan, Kishan Ganj, Delhi-7

Mob.: 9811256040

**INDOOR AED CABINET,
MODEL: Without Lock With Alarm**

Kay Ess AED cabinets with Alarm are made of sturdy heavy-gauge steel. AED Defibrillator Cabinet With Alarm can proudly display your AED for all to see and have the peace of mind knowing that it is protected from dirt, dust. In addition, the AED wall cabinets help the rescuers locate the AED device and provide an easy way to access it in case of an emergency.

Dimension



MANUFACTURED BY

KAY ESS TRADERS

Email: ravin.kayestraders@gmail.com

382-A, Gali No.4, Bagh Kare Khan, Kishan Ganj, Delhi-7

Mob.: 9811256040

INDOOR AED CABINET, MODEL: Without Lock With Alarm

Wireless Alarm

Clear PC Glass

**High
Quality
Magnet**

**Hanging
Accessories**

Safety Belt

**Door
Handle**

Key Features:

- AED Defibrillator Cabinet with Alarm.
- High Quality Magnet for keep the door close.
- Powder coated steel body for strength.
- The loud alarm emits a sound when the door is opened.
- Clear PC glass window and red & black lettering to facilitate visual identification of the AED.
- Economical Design with Smooth Finish.
- Easy to Install – Four holes are punched in the back of this cabinet for a wall mount.
- Three Side AED Graphics.
- Accommodates most AED's on the market.
- Longest Battery backup; Batteries Included.
- Accessories: Screws & plastic sleeves for fitting.
- Safety Belt.
- Product Dimension in Inches L x W x H: (6.5x15.5x16.5)
- Packing Dimension in Inches L x W x H: (9.5x18.5x18.75)
- Colour: White



**Wireless Alarm
On Off Switch**

