

MANUFACTURED BY

**KAY ESS TRADERS**

Email: ravin.kayestraders@gmail.com

382-A, Gali No.4, Bagh Kare Khan, Kishan Ganj, Delhi-7

Mob.: 9811256040

**INDOOR AED CABINET,  
MODEL: Without Lock Without Alarm**

Kay Ess AED cabinets without Alarm are made of sturdy heavy-gauge steel. AED Defibrillator Cabinet Without Alarm can proudly display your AED for all to see and have the peace of mind knowing that it is protected from dirt, dust. In addition, the AED wall cabinets helps the rescuers locating the AED device and provide an easy way to access it in case of an emergency.

### *Dimension*



MANUFACTURED BY


**KAY ESS TRADERS**

Email: ravin.kayesstraders@gmail.com

382-A, Gali No.4, Bagh Kare Khan, Kishan Ganj, Delhi-7

Mob.: 9811256040

**INDOOR AED CABINET,  
MODEL: Without Lock Without Alarm**



**Key features:**

- AED Defibrillator Cabinet Without Lock Without Alarm.
- High Quality Magnet for keep the door close.
- Powder coated steel body for strength.
- ClearPC glass window and red & black lettering to facilitate visual identification of the AED.
- Economical Design with Smooth Finish.
- Easy to Install – Four holes are punched in the back of this cabinet for a wall mount.
- Three Side AED Graphics.
- Accommodates most AED's on the market.
- Accessories: Screws & plastic sleeves for fitting.
- Safety Belt.
- Product Dimension in Inches L x W x H: (6.5x15.5x16.5)
- Packing Dimension in Inches L x W x H: (9.5x18.5x18.75)
- Colour: White

